

Quick & Easy Installation

1P_{ACK}

Can Be Ordered As A Kit With A General Use Seat Cushion



Lumbar Support General Use Wheelchair Back Cushion

Item #	Product Dimensions	Product HCPCS	UOM
14906	16" (W) X 17" (H) X 2.5" (D)	E2611	1/ea
14889	18" (W) X 17" (H) X 2.5" (D)	E2611	1/ea
14920	20" (W) X 17" (H) X 2.5" (D)	E2611	1/ea
8033	22" (W) X 19" (H) X 2.5" (D)	E2612	1/ea

- Memory-foam back cushion provides support, comfort and a tactile curve to improve posture.
- Durable Masonite board prevents hammocking in sling-back wheelchair. Easily removed as per patient preference.
- Conforms to lumbar curve to support lower back.
- Sealed liner protects foam core from moisture.
- Non-slip vinyl back prevents the cushion from sliding and is waterproof.
- Adjustable Straps easily attach to wheelchair.
- Meets flammability standard TB117.

specifications

Top Fabric Material	Black PU (stretch fabric)
Bottom Fabric	PVC + Knitted Fabric
Warranty	18 Months Limited