

wenzelite





HCPCS E1236



Trotter Mobility Chair

ltem #	Product Size	Product Seat Dimensions	Product Back Height	Product Seat-to Footplate	Product Seat-to- Floor	Product Weight	Product Weight Capacity	Product Weight Capacity used in Transit	UOM
TR 1200	12"	11.75" (W) x 11.5" - 13.5" (D)	20"	8"- 22"	22"	19 lbs.	75 lbs.	75 lbs.	1/bx
TR 1400	14"	13.75" (W) x 13" - 15" (D)	23"	9"- 22"	22"	21 lbs.	100 lbs.	100 lbs.	1/bx
TR 1600	16"	15.5" (W) x 13.5" - 15.5" (D)	24"	10"- 22.5"	22"	24.5 lbs.	170 lbs.	130 lbs.	1/bx
TR 1800	18"	17.5" (W) × 15" - 17" (D)	25"	11"- 23"	22.5"	30 lbs.	250 lbs.	130 lbs.	1/bx

- Lightweight folding frame with closure strap that holds the chair shut when folded.
- Black removable and washable fabric.
- Crash tested for transit. The Trotter has passed a dynamic crash test done with a 130 lb. front-facing occupant tested at 30 mph.
- Adjustable seat angle 15° and 22°.
- Adjustable back angle 85°, 90° and 95°.

specifications

Shipping Carton Weight (TR 1200)	28 lbs.
Shipping Carton Weight (TR 1400)	30 lbs.
Shipping Carton Weight (TR 1600)	33.5 lbs.
Shipping Carton Weight (TR 1800)	39 lbs.
Warranty	Limited Lifetime

- 8" front casters and 10" rear wheels with toggle lock.
- Adjustable-seat depth.
- Curb-assist lever facilitates curbside navigation. (Figure A)
- Positioning pelvic belt and five-point harness with padded strap covers. The reversible padded strap covers increase comfort and one side of cover is imitation lamb's wool and the other is velour.
- Swing-away footrests are independently adjustable and removable.