# SPECIFICATIONS

additional specifications at zingstanders.com

Weight Limit	<b>Size 1-</b> 70 lbs. (32kg) <b>Size 2-</b> 154 lbs. (70kg)
Height Range	Size 1- up to 44" (112cm) Size 2- 40"-60" (102-152cm)
Transfer Height	<b>Size 1-</b> 32"(81cm) <b>Size 2-</b> 31" (79cm)
Hip to Foot Plate Range	Size 1- 7"-22" (18-56cm) Size 2- 15"-33" (38-84cm)
Upper Body Support Height Range	Size 1- 4" (10cm) Size 2- 6" (15cm)
Knee Pad Depth Range	Size 1- 5" (13cm) Size 2- 7" (18cm
Weight of Basic Unit	<b>Size 1-</b> 55 lbs (25kg) <b>Size 2-</b> 82 lbs (37kg)
Frame Footprint	Size 1- 20"x28" (51x71cm) Size 2- 26.75"x39.75" (68x101cm)
30° Independent Leg Abduction on each leg	Available on all models
Supine	0°-90° (not available on Vertical & Prone)
Prone	0°-20°-MPS, MPS TT; 0°-90°-Prone (not available on Supine, Supine TT, & Vertical)
HCPCS	E0641-Multi-Position Stander MPS, MPS TT
	E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone
Positions	Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)
Adjustable Tray	Available on all models
Swing-Away Tray	Available (not available on Prone and Vertical)
Support	Individual Positioning Components, Padded Strapping
Dual Control	Available (not available on Vertical)
Gas Spring Lockout	Available (not available on Vertical)
Transfer Aid	Flat-to-Load (not available on Vertical)
Tool Free Adjustment	Yes
Pow'r Up Lift Option	Size 2 Only

#### **Full Line of Standers**

- MPS TT (tilt table)
- Supine
- Prone





## **Hygienic Covers**

- Available in 5 fun patterns for the pelvic, upper body, lateral, head and tray areas
- · Prevents stander from getting soiled, especially during feeding
- · Easily removed and washable





**SEA CREATURES** 







**Frame Colors** 










Altimate Medical, Inc. 262 West First Street Morton, MN 56270 USA Tel: 800.342.8968 507.697.6393 Fax: 877.342.8968 507.697.6900 info@zingstanders.com zingstanders.com

Visit our website at zingstanders.com for complete specifications on all Zing products and options.









# The 7ing Advantage









#### **True Multi-Position Stander**

- The MPS and MPS TT allow supine, vertical and prone positioning without flipping foot plates or pads, moving the tray or transferring out and back in the stander
- · Smooth transition from a flat-to-load transfer position, supine, vertical and 20 degrees prone - all in one motion, thanks to the unique pivoting action of the Zing MPS and MPS TT
- Provides therapeutic opportunities not possible with other multi-position standers - use prone for head control and when fatigue sets in, simply move to a slight supine position to rest but remain in the standing position
- Frame will only go supine and upright with Prone Position Lockout engaged

#### **True Hip Pivot Leg Abduction**

- Pivot points are in line with the hips, promoting proper alignment in the hips, knees and feet
- · Feet maintain the angle with the rest of the legs no wedges needed and no additional adjustment to the knees or foot plates
- Up to 30 degrees of abduction per leg accomplished through a single knob
- Promotes musculoskeletal development for children during critical growth years
- Abduction can be adjusted and used throughout the entire positioning range

#### Plenty of Growth

- Size 1: Fits children from infancy to 44", and up to 70lbs
- Size 2: Fits children from 40"-60", and up to 154lbs
- · Can easily grow with the child and adapt to their changing needs
- Many supports, pads and straps are available in multiple sizes
- Comfortable transfer height of 32" makes transfers more manageable as children grow

#### Remarkably Modular

- Can be configured as a simple stander or fit a child with the most complex needs
- Designed to grow with the child options can be added as needed
- · Great for a multi-user facility or school

# SUPINE TO PRONE











**Multi-Use Trav** 

Supine Lying - Flat-to-Load

- The MPS, Supine, MPS TT and Supine TT feature a multi-use tray that is angle adjustable and swings away either direction
- Tray can remain attached during transfers
- More activities are possible by swinging tray away during standing therapy

#### **Unique Lift Mechanisms**

- · Dual controls are available allowing angle adjustment by foot pedal or hand control
- · Gas Spring Lockout prevents an inadvertent change of position, especially important when around other children
- The Pow'r Up Lift option on Size 2 brings the user to the standing position with the touch of a button

### **Accommodating Contractures**

- No special modifications are needed to accommodate children with contractures
- Using the multi-adjustable foot plates, multi-adjustable knee pads and the independent calf pads, it is possible to adjust for contractures in one or both sides

## **Easy Adjustment**

- Tool-free adjustment of all components
- Every adjustment can be made with a simple turn of a knob or ratchet handle







