

# Protekt® Aire 7000

## Lateral Rotation/Low Air Loss/Alternating Pressure and Pulsation Mattress System



### FEATURES AND BENEFITS:

- Lateral rotation offers three turning modes with adjustable angles and turning cycle times from 10-95 minutes – Left Turn, Right Turn and Full Turn.
- TrueCare™ Alternating technology enhances patient comfort by re-inflating the 1st set of alternating cells before the 2nd set of cells start to deflate and alternate.
- Innovative iLAL™ (independent Low Air Loss) effectively reduces the moisture and skin heat on the support surface by periodically providing a high volume of air loss relief.
- Selectable 10, 15, 20, 25, 30 minute alternating cycle times with LED display and static mode for easy patient transfer.
- Constant low pressure (CLP) mode puts the mattress in a low pressure static surface for patients who cannot tolerate a dynamic surface but still require a pressure redistribution surface.
- Pulsate mode with adjustable cycle time provides a periodic pressure changing static surface.
- 18 bladder cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- Auto-firm mode provides quick set up and ease of ingress and egress from bed.
- Audio/visual alarms: power failure, low pressure and service alarm with error code.
- Nylon quilted top cover is fluid resistant, breathable, vapor permeable and fire retardant.
- Two-stage fowler detector automatically stops turning at 20 degrees and increases the pressure output at 45 degrees to prevent “bottoming out.”
- Auto-lockout ensures patient safety by locking the panel to prevent accidental setting changes.
- Weight Capacity: 660 lbs.
- Warranty: 2 Year Non-Prorated.

80070

Protekt® Aire 7000 36"x80"x8"

80070-42

Protekt® Aire 7000 42"x80"x8"



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com