



Transformer Scooter FAQ

What are the dimensions unfolded? **38"x17"x 25"**

What are the dimensions folded? **17.5" x 17" x 25"**

Lowest Tiller height from floorboard? **17.5"**

Highest tiller height from floorboard? **25.25"**

What are the body construction materials? **Steel Frame, Plastic Trim Pieces**

What is the height from top of the floorboard to ground? **6"**

What is the seat thickness? **3"**

What is the seat made of? **Foam, Vinyl covering**

What is the seat depth? (front to back) **15"**

What is the seat width? **17"**

What are the tires made of? **Solid rubber**

What is the ground clearance (bottom of floorboard)? **3.25"**

Distance from front of seat to floorboard? **13"**

Distance from front of seat to tiller? **11.75" to 13.75"**

How long does it take the battery to charge? **5-6 hours**

Can 2 remotes be programmed at the same time? **No. Only one can be synced.**