

SIMPLY OUTSTANDING

Stanzme is designed for daily standing therapeutic needs of children from early intervention to age 14. Its simple structure and application is easy to use by parents who may not have enough experience with this type of equipment. It is simple in design and construction as compared to others and allows caregivers to position the child quickly and with ease.



SUPINE STANDER

ADDRESS

16 Signal Hill Rd
Brookfield, CT 06804

TELEPHONE

+(888) 714-8405

EMAIL

sales@stanzme.com

WEB

www.stanzme.com



Stanzme Supine stander

The Stanzme Supine Stander is a versatile, easy to adjust option for children who need support in weight bearing and for whom prone standing is not an option. The child is always positioned at the lowest possible setting regardless of their size, and the stander “grows” with the child. A range of angle adjustments allows for a horizontal position so that the child can be transferred easily. The stander can then be raised to the desired angle.

Features

Basic stander includes:

- Lateral trunk and hip pads
- Padded back board
- Trunk and hip straps
- Lateral knee supports w/- straps
- Shoe holder w/- ankle and toe straps
- Tray with adjustable angle and removable bracket
- Locking casters
- Anterior chest harness
- Headrest with lateral adjustable pads



Supine Stander - Small



Supine Stander - Large



Specifications

Item	SSS-6000S	SSS-6000L
User Age/Weight	1 - 7 years/ 70lbs	7 - 14 years/ 130lbs
User Height (approximate)	33" - 55" / 84cm - 140cm	53" - 65" / 135cm - 165cm
Footrest dimensions	22"W x 12"D / 56cm x 32cm	25"W x 16"D / 64cm x 40cm
Base-width x length (fully extended)	21.5" x 41.5" / 55cm x 105cm	29.5" x 45.5" / 75cm x 116cm
Footrest to top of supine back	33.5" x 55" / 85cm x 140cm	53" x 59" / 135cm x 150cm
Tray Size	17" x 24" / 44cm x 61cm	19" x 28" / 48cm x 70cm
Range to Supine	Supine to 90°	Supine to 90°

Accessories



Two pairs of trunk and hip support offer alignment and comfort while improving and stabilizing balance and posture.



The tray is both height and angle adjustable. To adjust the height and angle, loosen the respective knobs located at the rear spine of the stander. The tray is removable for transfers.



Flat headrest with laterals is adjustable to provide midline head support.



Shoe holders with ankle and toe straps allow feet to stay in a secure and stable position.



Two 2" velcro straps are provided for trunk and hip supports to securely position the child in place.



Locking casters prevent roll and swivel during transfers.



Lateral knee supports ensure the knees are properly supported at the correct height for each child.

Testimonial



As an ATP with over 25 years of experience in custom rehab equipment and standing equipment solutions, Stanzme is my choice of therapeutic equipment as it is simple and easy to use. A child is transferred on it with ease by the caregiver, and secured with prompts from arms, to chest, to hips, to knees, to feet that are easily adjustable.

I like the product because of its ease of use and simplicity, and the child benefits from daily standing therapy. I've no hesitation in recommending Stanzme stander as it can be used efficiently and as often as possible without any technical knowledge or assistance from therapists on a regular basis. It is simple, affordable and insurance-friendly product designed to provide maximum comfort with minimum hassle.

- Joe Q, ATP

Special discount price for Dealers and Special Education Schools. Call us: +(888) 714-8405

ADDRESS

16 Signal Hill Rd
Brookfield, CT 06804

TELEPHONE

+(888) 714-8405

EMAIL

sales@stanzme.com

WEB

www.stanzme.com