



**FUTURE  
MOBILITY**



**CAPELLA**  
BY FUTURE MOBILITY

# TILT WHEELCHAIR

# THE FUTURE OF MOBILITY – WORKING FOR YOU

"Assembled, Owned and made in North America." In 2000 Future Mobility began designing, manufacturing, and then distributing wheelchairs for the North American Home Healthcare market. Our Founder's mission is to help the world by designing high quality, durable mobility products for the North American marketplace.



ORION II

ORION II 500

ORION III

SUPERKIDS



## DYNAMIC TILT

OPTIONS: 20 degrees of tilt (optional 45 degrees of tilt with Capella 45), 21 degrees of recline (optional), growable frame – Seat width can increase from 16" to 20" or 20" to 24". The seat's depth can increase by as much as 1" without the need to change frame members.



### COLOR OPTIONS



## LOWEST SEAT-TO-FLOOR-HEIGHT ON THE MARKET




<b>40 LBS. WEIGHT</b>	<b>20° &amp; 45° TILT DEGREE</b>	<b>350 LBS. CAPACITY</b>
---------------------------	--------------------------------------	------------------------------

This Dynamic Tilt is designed with the lowest Seat-To-Floor Height on the market 11.75". It's Ideal for foot propellers



Scan me






## 11.75" SEAT-TO-FLOOR HEIGHT

## QUICK RELEASE

Allows for quicker and easier changing or removing of the wheelchair tires. This option is of value to those who use the tilt chair for transport and want to reduce weight of the chair when loading it into a vehicle.

## ADJUSTABLE CASTER HEIGHT

There are different types of caster height-to floor options available. Select the right height according to the requirement.




## ABILITY TO GROW

Frame grows as required without replacing any parts up to 4" of width growth. Seat depth is growable up to 2" with each seat rail at no extra cost.

## GAS CYLINDER

Single cylinder linkage mounts onto dual slide rails for ultra rigid frame and simplified design. It can be found located underneath the chair, allows for a gentle, slow tilt for the client lockable at infinite angles






# SPECIFICATIONS:

**CAPELLA**  
BY FUTURE MOBILITY  
**DYNAMIC TILT**




**SEAT TO FLOOR**  
11 3/4" - 19"  
(13" - 19" ON CAPELLA 45°)



**REAR WHEEL**  
12" - 16", 20", 22" OR 24"  
COMPOSITE URETHANE OR  
SPOKED WHEELS



**SEAT WIDTH**  
16" - 24"



**SEAT DEPTH**  
15" - 20"



**OVERALL WIDTH**  
11 3/4" - 19"  
(13" - 19" ON CAPELLA 45)




**CASTER SIZE**  
4" (ONLY ON CAPELLA)  
5", 6", 7", 8" URETHANE



**FRONT RIGGING**  
SWINGAWAY  
DUAL SWINGAWAY  
OR DUAL ELEVATING LEGRESTS



**PRODUCT WEIGHT**  
66 LBS (30KG)




**ADJUSTABLE  
BACK ANGLE**  
-7°, 0°, 7°, 14°, 21°



**BACK STYLE**  
ORON BACK



**ARM STYLE**  
ADJ. HEIGHT T-STYLE  
ADJ. HEIGHT PIN RELEASE T-STYLE  
ADJ. HEIGHT FLIP BACK



**FRAME MATERIAL**  
ALUMINUM



**WEIGHT CAPACITY**  
CAPELLA 20° - 350 LBS (159 KG)  
CAPELLA 45° - 250 LBS (113 KG)

**ADP**

**CODE**  
WM0000559

**HCPCS**

**CODE**  
E1161

## WARRANTY:

The Capella, Capella45 is warranted to be free from defects in materials and workmanship for a period of one (1) years and the Orion II 500 for period of one (1) year from the date of purchase. The frame has a lifetime warranty and the upholstery for 90 days, upon normal usage by the original purchaser. Power tilt batteries are warranted for six (6) months while the control box, battery charger, actuator and hand control are warranted for one (1) year. If within this warranty period the product shall be proven to be defective, such product shall be repaired or replaced, at Future Mobility Healthcare Inc. discretion.

**USA**



11126 Greenstone Ave, Santa Fe Springs, CA 90670



716-783-9130



[www.futuremobility.com](http://www.futuremobility.com)



**FUTURE MOBILITY**

**CANADA**

6750 Professional Court, Mississauga, Ontario, L4V 1X6



1-888-737-4011



[www.futuremobility.ca](http://www.futuremobility.ca)

