







CH4808-L-BLUE 870 Regulator

For use with: M2 through E Cylinders

• Ergonomic Handle



CH3115-L 870 Regulator

For use with: M2 through E Cylinders

Brass Core

CHAD® Regulators

Item #	Product Description	Settings	Outlet	Connection	Liter Flow Increments	UOM
CH4808-L-BLUE	CHAD® Click Style Blue Regulator	0-8 lpm	Barb	CGA 870	0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8	1/ea
CH4815-L-BLUE	CHAD® Click Style Blue Regulator	0-15 lpm	Barb	CGA 870	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	1/ea
CH3115-L	CHAD® Brass Core Regulator	0-15 lpm	Barb	CGA 870	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	1/ea
CH3215-L	CHAD® Brass Core Regulator	0-15 lpm	Barb	CGA 540	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	1/ea
CH3104-L	CHAD® Brass Core Pediatric Regulator	0-4 lpm	Barb	CGA 870	1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1.5, 2, 3, 4	1/ea

CHAD® Click Style and Brass Core Regulators

- 20-micron, sintered-bronze inlet filter.
- FDA-recommended viton and brass seal washers.
- All flows tested to accuracy within +/- 15% tolerances.
- Meets ASTM-G175-03 promoted ignition tests.
- Inlet pressure: 200 3,000 PSI.
- Outlet pressure: 50 PSI.
- 100% QC tested in the USA.

specifications	
Warranty	Limited 3 Year

accessories		UOM
SEAL3000	Brass & Viton Triangle Seal Washer for Oxygen	10/pkg
	Regulators & Conservers	