

Deluxe Clinical Care Recliner with Drop Arms

The Deluxe Clinical Care Recliner sets a new standard of comfort for the patient and ease of use for the caregiver.

- Innovative, easy to actuate drop-arm feature. Self-storing drop-arms are simple for the caregiver to activate and enables a safer, easier lateral transfer for the resident and caregiver.
- Available With Heat and Massage Option.
- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Easy to actuate by caregiver. Side table surface is smooth and features tapered edges for easy cleaning.
- The Series FR566G and FR566GH Deluxe Clinical Care Recliners feature three positions for patient comfort. The fourth position is a Trendelenburg (or Shock Position) that must be activated by a caregiver. The Trendelenburg Position cannot be activated by a patient while in the recliner.
 - The first position is the **Upright Position**. This position is ideal for patient transport and can be used to initiate patient treatment / therapy.
 - The second position is the **First Recline Position** and is also known as the “TV” position. This position is ideal for patient treatment, rest and reading.
 - The third position is the **Second Recline Position** and is also known as the “Full Recline” position. This position is ideal for patient treatment, rest and reading.
- The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of comfort. Built-in arm pads in the back give improved support for the arms.
- Wall saver armrest caps increase durability of the armrests and protect walls from marring.
- Higher floor to back height enables the caregiver to raise the patient up with less effort, reducing the chance for a back injury.
- Extra deep Trendelenburg Position is easily actuated from either side of recliner by a conveniently located foot pedal. Steeper angle helps improve patient stabilization.
- Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with easy to locate black locking tab lock wheel and swivel with single motion. Directional wheel lock with high visibility green tab allows recliner to track straighter with less effort during transport.
- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Easy to actuate by caregiver. Side table surface is smooth and features tapered edges for easy cleaning.
- Extra heavy duty recliner mechanism.
- 350 lb maximum weight capacity, EVENLY DISTRIBUTED.



Specifications

Overall Height.....	51.5"
Overall Width.....	33"
Overall Depth.....	38"
Width Between Arms.....	23"
Arm Height From Floor.....	28"
Arm Height From Seat.....	8"
Seat Height.....	22"
Seat Width.....	20.5"
Seat Depth.....	19"
Depth Fully Reclined	
In Trendelenburg.....	76"
Width With Both Side Tables	
in Upright Position.....	53"
Shipping Weight.....	135 lb
Max. Weight Cap.....	350 lb



Multiple Positions for Enhanced Patient Comfort



Colors may vary slightly from actual vinyl

Standard

with Heat & Massage

FR566DG401	FR566DGH401	Vintage Gold
FR566DG409	FR566DGH409	Warm Taupe
FR566DG427	FR566DGH427	Blue Ridge
FR566DG432	FR566DGH432	Imperial Blue
FR566DG454	FR566DGH454	Royal Blue
FR566DG6705	FR566DGH6705	Dolce Sand
FR566DG6714	FR566DGH6714	Dolce Burgundy
FR566DG6717	FR566DGH6717	Dolce Moss
FR566DG6724	FR566DGH6724	Dolce Jet
FR566DG6725	FR566DGH6725	Dolce Indigo
FR566DG6727	FR566DGH6727	Dolce Meteor
FR566DG8504	FR566DGH8504	Sand
FR566DG8516	FR566DGH8516	Chestnut
FR566DG8533	FR566DGH8533	Port
FR566DG8561	FR566DGH8561	Mulberry
FR566DG8567	FR566DGH8567	Dove
FR566DG8584	FR566DGH8584	Taupe
FR566DG8602	FR566DGH8602	Oatmeal
FR566DG8604	FR566DGH8604	Concrete
FR566DG8611	FR566DGH8611	Terra Cotta

Standard

with Heat & Massage

FR566DG8615	FR566DGH8615	Berry
FR566DG8624	FR566DGH8624	Granny Smith
FR566DG8626	FR566DGH8626	Eucalyptus
FR566DG8628	FR566DGH8628	Sea
FR566DG8801	FR566DGH8801	Marine Blue
FR566DG8805	FR566DGH8805	Jet
FR566DG8808	FR566DGH8808	Champagne
FR566DG8809	FR566DGH8809	Taupe
FR566DG8818	FR566DGH8818	Sage
FR566DG9201	FR566DGH9201	Cobblestone
FR566DG9202	FR566DGH9202	Periwinkle
FR566DG9205	FR566DGH9205	Gypsum
FR566DG9207	FR566DGH9207	Ice Blue
FR566DG9208	FR566DGH9208	Wineberry
FR566DG9210	FR566DGH9210	Cypress
FR566DG9212	FR566DGH9212	Blue Jade
FR566DG9214	FR566DGH9214	Steel Blue
FR566DG9715	FR566DGH9715	Eucalyptus
FR566DG9722	FR566DGH9722	Desert Sage
FR566DG9764	FR566DGH9764	Steel Grey