


UNIVERSAL THERAPY SYSTEM

Universal Therapy System is a low to medium risk system for patients who are at risk of developing pressure ulcers according to your sound clinical judgment. UTS can also be used for patients who have an existing stage 1, 2, 3 and 4 pressure ulcers, in conjunction with your hospital's policy on pressure area management.



INTEGRATED CONTROL UNIT	
Dimension:	(L)12.1"x (W)7.5"x (H)5.3"
Weight:	7.1 lbs
Rating Voltage:	110-120 volts
Max Current:	0.1A max.
Rating Frequency:	60Hz
Cycle Time:	10 mins
	Static function
	2 in 1 alternation
Min/Max Pressure:	12~47 mmHg +/- 5 mmHg
Flow-rate:	> 5 l/min
Protection Type:	Class I, Type BF
Low Pressure Alarm:	Visual / Audio
Mode Selection:	Static/Alternation Function

MATTRESS REPLACEMENT	
Model No:	Universal Therapy System
Dimensions:	(L)80"x (W)36"x (H)7"
Weight:	31 lbs
Weight Capacity:	350 lbs
Cell Material:	Nylon with PU coating
	Foam in air cell design for protection during power outage
Cell Number:	21 cells
Cover Material:	Stretch Nylon with PU coating
Base Material:	Polyester with PVC coating
	CPR Quick Release System
HCPSC Code	E0277

PRIUS HEALTHCARE USA
Creating innovation and difference

T (813) 854-5464
E info@priushc.com
W priushcusa.com

160 SCARLET BLVD
OLDSMAR, FL 34677

UNIVERSAL THERAPY SYSTEM


Provides you, the caregiver, with peace of mind, and your patient with world class care and comfort



PRIUS HEALTHCARE USA
Creating innovation and difference

priushcusa.com


Perfect for patients and caregivers



PRESSURE MAPPING

UTS provides constant low pressure by redistributing the patient's pressure over a greater area than a standard mattress.

UNIVERSAL THERAPY SYSTEM



- The Universal Therapy System has an integrated pump that is automatic weight sensing
- Alternating pressure and continuous minimized pressure therapy modes to prevent tissue breakdown and encourage healing
- Foam in air cell design provides safe and sound protection during a power outage
- Removable foam layer that can be placed under top cover for extra comfort or stored underneath the air cells
- Hand pendant for patient control of the mattress. Ideal for patients living alone or with limited mobility.
- The UTS replacement mattress can accommodate patient weights up to 350 lbs

Universal Therapy System is your partner in protection and treatment of pressure ulcers



Low friction/low shear 2-way stretch top cover



CPR quick release system removes the air within 20 seconds

Hand pendant for patient control



Welded vinyl bladder containment straps



Protective foam pad covers the integrated pump

Hand Pendant Controls

- Power On/Standby Selection
- Static/Alternating Mode
- Comfort Setting
- Automatic Weight Sensing
- Function Status Display
- Max Inflate
- Alarm Display
- Lock/Unlock & Alarm Mute



Hard-plummed air cells to prevent leaks

