

BARIATRIC SHOWER COMMODE CHAIR



Gendron's new shower commode chair features tensioning back support straps to better position the individual. The back upholstery allows water to pass through and is mildew and bacteria resistant. Arm rests are removable for lateral transfer. The commode pail is removable from the rear of the chair. Two locking rear casters. The frame is constructed of rust resistant powder coated stainless steel. Full width push bar at the top of the back.



Model	Seat Width	Seat Depth	Seat to Floor Ht.	Arm Style	Overall Width	Footrests
5226	26"	20"	21.7"	removable	35.3"	removable
5228	28"	20"	21.7"	removable	35.3"	removable

GENDRON 
Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
 Phone: 800.537.2521 www.gendroninc.com